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Information for optimization of the work result:

- Best results are achieved by working with slight pressure as possible and with brushing movements. Please comply with the working pressure of a few pond for all FG and RA instruments up to diameter 023. Thus, work with as low pressure as possible. RA diamond burs with diameter 025 and more are resilient to up to 500 pond of working pressure.
- After grinding with coarse and extra-coarse grit, finishing of the surface is required.

Processing Information

Storing

Please store all originally packaged tools in a clean and dry place and at room temperature. All rotary instruments must be disinfected, cleaned and - if applicable - sterilized prior to the first and after each subsequent use.

Cleaning

- For pre-cleaning of diamond instruments, we recommend the use of DIACLEAN cleaning stone (Order no. 500).
- Cleaning may be done in a cleaning machine, as long as appropriate clamping implements and cleaning agents are used that are suitable for rotating instruments.
- Diamond instruments are suitable for ultrasonic baths as well, if appropriate clamping implements are used, which prevent touching of the heads.
- For disinfection, disinfection baths can be used. Please use disinfectants with corrosion protection only and observe the contact times prescribed by the manufacturer.

Sterilization

- Polishers, diamond and carbide instruments can be vapour-sterilized in an autoclave at up to 135°C. We recommend the following sterilization parameters: 3x fractionated pre-vacuum, sterilization temperature 134°C, holding time 5 min and drying time 10 min*.
- Please make sure, e.g. by using suitable clamping implements, that the blades or diamond heads are not touched.
- Polishers must not be heated above 135°C, and the cleaning and disinfectant must be explicitly suitable for them.
- Carbide instruments must not be heated above 180°C nor exposed to hydrogen peroxide or strong acids.
- Heat or hot air sterilization and thermal disinfection are not suitable for our instruments!

Please follow precisely the instructions provided by the manufacturer!

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